



ANESTHESIA INSTRUMENTS

✉ info@sacmed.bz
🌐 www.sacmed.bz
☎ +92 52 3256007
📞 +92 321 612 7007
📍 Noul Mour, Opp.
Anwer Khawaja Roras Road,
Sialkot - 51310, Pakistan.


- Low Profile
- These are single-patient-use alternatives to our reusable blades. Along with high-quality handles & fiber optic technology, it gives the same feeling as reusable ones.
- Single Use blades help eliminate the time, reprocessing steps and cost related to disinfecting reusable blades.
- Reduced risk of cross-contamination
- Stainless steel body
- Non-Sterile, clinically clean product as per FDA guidelines
- Fiber Optic Light



Product Code	Blades	Fig.	Working Length
MC0-504-630	Neonate	0	57 mm
MC1-504-631	Infant	1	72 mm
MC2-504-632	Child	2	92 mm
MC3-504-633	Medium Adult	3	113 mm
MC4-504-634	Large Adult	4	140mm

Non-Sterile | Unit of Sale 20 Pcs / Box

- Low Profile
- These are single-patient-use alternatives to our reusable blades. Along with high-quality handles & fiber optic technology, it gives the same feeling as reusable ones.
- Single Use blades help eliminate the time, reprocessing steps and cost related to disinfecting reusable blades.
- Reduced risk of cross-contamination
- Stainless steel body
- Non-Sterile, clinically clean product as per FDA guidelines
- Fiber Optic Light



Product Code	Blades	Fig.	Working Length
ML00-504-729	Premature	00	47 mm
ML0-504-730	Neonate	0	57 mm
ML1-504-731	Infant	1	77 mm
ML2-504-732	Child	2	132 mm
ML3-504-733	Medium Adult	3	172 mm
ML4-504-734	Large Adult	4	182 mm

Non-Sterile | Unit of Sale 20 Pcs / Box

- Single-use blade with attached handle (Handle and blade fixed together cannot be separated)
- Low profile
- Pre-loaded with batteries
- Surgical grade stainless steel blade
- High-grade, reinforced aluminum handle with an ergonomic grip
- Fiber optic light provides optimal light output
- Color-coded sizes for easy identification
- Non-Sterile, clinically clean product as per FDA guidelines



Product Code	For	Fig.	Working Length	Color
MCCO-504-6607	Neonate	0	54 mm	
MCC1-504-6617	Infant	1	69 mm	
MCC2-504-6627	Child	2	89 mm	
MCC3-504-6637	Medium Adult	3	110 mm	
MCC4-504-6637	Large Adult	4	137 mm	

Non-Sterile | Unit of Sale 10 Pcs / Box

- Single-use blade with attached handle (Handle and blade fixed together cannot be separated)
- Low profile
- Pre-loaded with batteries
- Surgical grade stainless steel blade
- High-grade, reinforced aluminum handle with an ergonomic grip
- Fiber optic light provides optimal light output
- Color-coded sizes for easy identification
- Non-Sterile, clinically clean product as per FDA guidelines



Product Code	For	Fig.	Working Length	Color
MLC00-504-7697	Premature	00	44 mm	
MLC0-504-7607	Neonate	0	54 mm	
MLC1-504-7617	Infant	1	74 mm	
MLC2-504-7627	Child	2	129 mm	
MLC3-504-7638	Medium Adult	3	169 mm	
MLC4-504-7648	Large Adult	4	179 mm	

Non-Sterile | Unit of Sale 10 Pcs / Box

- Single Use Laryngoscope Handles
- Material: Aluminium
- Single Patient use
- Individually Packed
- Clean, Non Sterile
- Clinically clean product as per FDA guidelines

Product Code

Size

SH-604-5121 Small / Penlight
MH-604-5120 Medium



Medium

Small / Penlight

Non-Sterile | Unit of Sale 20 Pcs / Box

- Design for use in 6.0mm and larger size endotracheal tube intubation placement.
- Easy to use, learn, and teach
- Clinically clean, packed in clean room and ETO sterilized
- Provides maneuverability for placement of an endotracheal tube
- Easy to high-level disinfect sterilized by low temperature ethylene oxide gas.
- Load the endotracheal tube onto the stylet, ensure that the distal end of stylet does not extend beyond the distal end of the ET tube.
- Position the ET tube at the opening of the vocal cords. Do not advance the stylet past the vocal cords.
- When the ET tube is engaged in the glottic openings, extract the stylet 5 cm by use of the thumb pushing up on the thumb tab of the stylet. The partial removal of the stylet from ET tube softens the tip, allowing the tube to pass through the vocal cords.
- Position the ET tube as per standard practice.
- Remove the stylet from ET tube.

Product Code

Ball Size

Working Length

Overall Length

SM-DPS125 5 mm Diameter 266 mm (10.47") 360 mm (14.17")



Sterile | Unit of Sale 25 Pcs / Box

Reusable / Autoclaveable


- Surgical grade stainless steel
- Smooth Edges
- Rounded tip & high quality surface finish
- Light guide design allow easier cleaning and extends blade life
- Provides with 4mm fiber optic for light pass gives true tissue color
- Sand finish

Product Code

605-111

Blade	Fig	Working Length
605-111-1	1	70 mm
605-111-2	2	91 mm
605-111-3	3	110 mm
605-111-4	4	131 mm
605-111-5	5	155 mm

Unit of Sale: 1 Set



Reusable / Autoclavable

- Surgical grade stainless steel
- Smooth Edges
- Rounded tip & high quality surface finish
- Large fiber bundles offer more light no black marks
- Fully autoclavable and easier to clean
- Sand finish

Product Code **605-210**

Blade	Fig.	Working Length
605-210-1	1	72 mm
605-210-2	2	92 mm
605-210-3	3	112 mm
605-210-4	4	135 mm

Unit of Sale: 1 Set



MILLER FIX TUBE HBT (FIBER OPTIC) COLD LIGHT GUIDE

Reusable / Autoclavable

- Surgical grade stainless steel
- Smooth edges
- Rounded tip & high quality surface finish
- Provides with 5 mm fiber optic offers no black marks
- Fully autoclavable and easier to clean
- Sand finish

Product Code **605-211**

Blade	Fig.	Working Length
605-211-0	0	57 mm
605-211-1	1	82 mm
605-211-2	2	132 mm
605-211-3	3	172 mm
605-211-4	4	181 mm

Unit of Sale: 1 Set



Reusable / Autoclaveable

- Surgical grade stainless steel
- Smooth Edges
- Rounded tip & high quality surface finish
- Provides with 5mm fiber optic for light pass gives true tissue color
- Fully autoclavable and easier to clean
- Sand finish



Product Code **605-214**

Blade	Fig.	Working Length
605-214-1	1	70 mm
605-214-2	2	90 mm
605-214-3	3	110 mm
605-214-4	4	131 mm

Unit of Sale: 1 Set



INTEGRATED LARYNGOSCOPE (FIBER OPTIC)

Reusable / Autoclaveable

- Surgical grade stainless steel
- Smooth Edges
- Rounded tip & high quality surface finish
- Provides with 5 mm fiber optic for light pass gives true tissue color
- Fully autoclavable and easier to clean
- Sand finish



Product Code **605-215**

Blade	Fig.	Working Length
605-215-0	0	57 mm
605-215-1	1	82 mm
605-215-2	2	132 mm
605-215-3	3	172 mm
605-215-4	4	181 mm

Unit of Sale: 1 Set



Reusable / Autoclaveable

- Surgical grade stainless steel
- Smooth Edges
- Rounded tip & high quality surface finish
- Spring ball bearings gives secure & firm support supplied with LED or Halogen bulb
- Sand finish

Product Code

605-216

Blade	Fig.	Working Length
605-216-1	1	72 mm
605-216-2	2	92 mm
605-216-3	3	112 mm
605-216-4	4	131 mm
605-216-5	5	156 mm

Unit of Sale: 1 Set



MILLER CONVENTIONAL

Reusable / Autoclaveable

- Surgical grade stainless steel
- Smooth Edges
- Rounded tip & high quality surface finish
- Spring ball bearings gives secure & firm support supplied with LED or Halogen bulb
- Sand finish

Product Code

605-217

Blade	Fig.	Working Length
605-217-0	0	58 mm
605-217-1	1	83 mm
605-217-2	2	113 mm
605-217-3	3	133 mm
605-217-4	4	182 mm

Unit of Sale: 1 Set

